

Untreated test tubes sterile

for example for the storage and transport of swab samples



Untreated test tubes 12 ml, sterile

Polystyrene, with pressed-on (overlapping) stopper
Ø 16 mm, L = 95 mm

Item No.	Description
050401	N 95 S sterile
050402	N 95 Sn sterile

Packing unit: 100 pcs. in bag, 1,000 pcs. in box




Untreated test tubes 12 ml, sterile

Polypropylene, with pressed-on overlapping stopper
Ø 16 mm, L = 95 mm

Item No.	Description
050502	N 95 S PPn sterile

Packing unit: 100 pcs. in bag, 1,000 pcs. in box



Untreated test tubes 13 ml, sterile

Polypropylene, with screwed-on cap
Ø 16.5 mm, L = 101 mm

Item No.	Description
0961 1150	N 101/16,5 S sterile

Packing unit: 100 pcs. in bag, 1,000 pcs. in box